


### Fully Sealed Edge (pouch pad)

Poly film top layer Non-woven PET bottom




### Multiple Structural Options

- All Paper
- SAP (Super absorbent Polymer)
- Wide range of absorbencies



### Perfect for any Tray Application

- Variety of sizes available
- Colors to match any tray
- Can be used in foam or solid sheet trays



**Pactiv** Absorbent Pads and Case Liners are designed to meet a wide variety of applications. From Protein to Agricultural Processors and Grocer to Medical applications, **Pactiv** is a trusted supplier for your Absorbent Pad and Case Liner Needs.

- Stock and customized sizes available
- Higher moisture retention than conventional pulp/paper pads
- Lighter, thinner pads
- Available in a variety of sizes, absorbencies and colors



# Stock Absorbent Pad and Case Liner Offerings

## All Paper

### Edge Sealed Absorbent Tray Pads

- Economical
- Absorbencies from 25-300g
- Standard white or custom colors

Description	Size (in.)	Absorbency (g)	Pads/Case
No Wet 25	3" x 7"	25	3,000
No Wet 40	4" x 7"	40	2,000
No Wet 50	4" x 7"	50	2,000

## SAP

### Edge Sealed Super Absorbent Polymer Tray Pads

- Retain 35-50% more under load than all paper pads
- Fresh, dry customer package
- Clean retail cases
- Standard white or custom colors



Description	Size (in.)	Absorbency (g)	Pads/Case
No Wet 25S	3" x 6"	25+	3,000
No Wet 40S	4.5" x 6"	40+	2,000
No Wet 50S	4.5" x 6"	50+	2,000

## Absorbent Case Liners

### Paper Ply Absorbent Core

- Absorbencies from 15g to 1000g
- CO2 generation option
- White, Black, White/Black color options
- Poly 1 or 2 side in singles sheets cut to size or roll form

Description	Size (in.)	Absorbency (g)	Pads/Case
White Poly 1 Side	30" X 50"	4 ply	50 Sheets
White Poly 2 Side	30" X 250"	4 ply	1 Roll
Black/White Poly	30" X 250"	4 ply	1 Roll



\* Tray colors, dimensions, availability may vary by region.

**Contact your Pactiv Sales Representative for more information.**

These guidelines are supplied to assist you in determining the proper use of Pactiv products. They are based upon testing and publishing guidelines and are reliable in most applications. However, because every food supplier's recipes, ingredients, processes and supply chain is unique, these guidelines are not a substitute for product testing. Confirmation of product acceptability under your specific conditions of use must be done by you.